



ROTIRAJUĆI DENTALNI INSTRUMENTI



**NAJVEĆI ASORTIMAN NA JEDNOM MESTU!**

[www.rotodent.rs](http://www.rotodent.rs)

## O NAMA

Roto Dent se bavi prometom rotirajućih instrumenata u stomatologiji i zubnoj tehnici.

- ▶ Od osnivanja zastupamo Austrijsko/Švajcarsku firmu **G&Z Instrumente GmbH** sa sedištem u gradiću Lustenau u Austriji. Sopstvena proizvodnja, kao i manualna kontrola kvaliteta izrade svakog proizvoda, dokaz su zašto se ova firma nalazi u samom svetskom vrhu proizvođača iz ove oblasti.
- ▶ Takođe, zastupamo i dobro poznatu Kanadsku firmu **TRI HAWK**, najpoznatijeg i najboljeg proizvođača karbidnih borera za sečenje keramike i metala.
- ▶ Početkom 2020. smo registrovali proizvod **PLASTOR** - privremeni ispun za zatvaranje kaviteta vrhunskog kvaliteta, **italijanskog proizvođača GHIMAS s.a.**

Nadamo se da ćemo kvalitetom i brzinom naših usluga uspeti da zadobijemo i Vaše poverenje.


S poštovanjem,

**Željko Uskoković**

# TEHNIČKI DETALJI

## Finoća

<b>UF</b>	ultra fino	15 $\mu$	<b>N</b>	srednje	107-126 $\mu$
<b>XF</b>	vrlo fino	20 $\mu$	<b>C</b>	grubo	151 $\mu$
<b>F</b>	fino	40 $\mu$	<b>SG</b>	veoma grubo	181 $\mu$

 FG Preporučeni broj obrtaja u minuti

$\phi$ in 1/10 mm	preporučeni (min-1)	maksimalni (min-1)
008 ° 010	75.000 ° 150.000	450.000
012 ° 014	60.000 ° 110.000	450.000
016 ° 018	45.000 ° 88.000	450.000
021 ° 023	40.000 ° 75.000	300.000
025 ° 027	30.000 ° 65.000	160.000
029 ° 031	25.000 ° 56.000	140.000
033 ° 040	22.000 ° 45.000	120.000
042 ° 050	20.000 ° 37.000	95.000

## Oblik tela borera

### FG

Za turbinu

$\phi$  1,6mm


16,5mm



19 - 24mm



21mm



25mm




### RA

Za kolenjak

$\phi$  2,35mm

22mm



26mm




34mm



# VIŠESLOJNI DIJAMANTSKI BORERI

✓ **Dijamantski boreri** mogu se koristiti za različite stomatološke restaurativne postupke, uključujući, ali ne ograničavajući se na **uklanjanje kariesa, starih stomatoloških restauracija i drugih tretmana**

✓ **Za optimalne performanse** uvek prilagodite model i finoću borera sa kojim radite



# VIŠESLOJNI DIJAMANTSKI BORERI




*Dijamantske  
čestice različitih  
veličina i finoća*




## ZLATNI BORERI

- ✓ Izuzetna trajnost
- ✓ Brza obrada i preciznost



## ZD2 BORERI

- ✓ Nova serija ZD2 borera, omogućava rad sa veoma tvrdim materijalima kao što je Cirkon oksid ( $ZnO_2$ )
- ✓ Fino doterivanje/brušenje gotovih krunica ili uklanjanje starih






## SPIRALNI BORERI

- ✓ Posebno konstruisan oblik tela borera
- ✓ Brži rad u poređenju sa standardnim borerima
- ✓ Bolje hlađenje i odnošenje viška materijala

# KARBIDNI BORERI

- ✓ Karbidni boreri su **nezamenjivi** u stomatološkoj praksi
- ✓ Omogućavaju **precizno i lako uklanjanje** kariesa
- ✓ Napravljeni su od čvrstog materijala koji omogućava **oštar rez** i dugotrajnu upotrebu
- ✓ Dostupni su u različitim oblicima - **okrugli, cilindrični ili konusni**



**Preciznost i brzina** smanjuju trajanje tretmana

**Izdržljivost** ih čini idealnim alatom za dugotrajnu upotrebu


Dostupni su u različitim oblicima i veličinama

**Manja vibracija** smanjuje osećaj nelagodnosti za pacijente

**Lako se čiste** i sterilizuju

# HIRURŠKI KARBIDNI BORERI

- ✓ Omogućavaju **precizno uklanjanje tkiva ili koštane strukture** pri izvođenju različitih hirurških tretmana
- ✓ **Veoma kvalitetni instrumenti** izrađeni od posebnog karbidnog materijala koji garantuje dugotrajnost, oštrinu i preciznost u radu



 Preporučeni broj obrtaja u minuti

# GUMICE ZA POLIRANJE

GUMICE ZA POLIRANJE



## GUMICE ZA AKRILAT



2234F

2236F

2238F

2230F

2130M

2134M

2136M

2030C

2034C

2036C

5.000 - 7.000 rpm  
max. 15.000 rpm

5.000 - 7.000 rpm  
max. 15.000 rpm

5.000 - 7.000 rpm  
max. 15.000 rpm

## KERAMIKA



TWIST  
GRUBA

TWIST  
PREDPOLIRANJE

TWIST  
POLIRANJE

TWIST  
GRUBA

TWIST  
PREDPOLIRANJE

TWIST  
POLIRANJE

OBRADA ZIRKON  
KERAMIKE

10.000 - 12.000 rpm  
max. 20.000 rpm

7.000 - 12.000 rpm  
max. 15.000 rpm

za poliranje i predpoliranje



# ČETKICE ZA POLIRANJE KERAMIKE I KOMPOZITA

sa silikon karbidom



WT1 WZ1 WS1

sa dijamantom



WTD1 WZD1 WSD1

# GUMICE ZA POLIRANJE

## UNIVERZALNE



UNI 1071 UNI 1092 UNI 1079 UNI 1077

7.000 - 12.000 rpm  
max. 20.000 rpm

## KOMPOZIT



TWIST MEDIUM TWIST SOFT ECC 5177M ECC 5188M ECC 5173M

7.000 - 12.000 rpm  
max. 20.000 rpm

7.000 - 10.000 rpm  
max. 20.000 rpm


## KOMPOZIT



DIACOMP ULTRA SA DIJAMANTOM

One set polishing

3 000 - 10 000 rpm  
Max 20 000 rpm



7.000 - 10.000 min<sup>-1</sup>  
max. 20.000 min<sup>-1</sup>

3.000 - 8.000 min<sup>-1</sup>  
max. 20.000 min<sup>-1</sup>

# KARBIDNE FREZE



79-60SG 79-70SG 79-40C 79-60C 251-60C 73-60N 78-60N 79-60N 251-60N 261-23N 138-16N



295-23N 77-60F 79-31F 79-45F 79-60F 257-23F 251-60F 295-23F 79-40TC 251-60TC 261-23TC



Najgrublja



Gruba



Standardna



Fina



Superfina

## PREPORUČENE BRZINE

010 - 023 5 000 - 40 000 rpm

025 - 035 15 000 - 30 000 rpm

HP

Handpiece

Ø 2,35mm

32mm

44,5mm

70mm

UM: unmounted

# DISKOVI



**ZD N 366**



**ZD F 366**



**ZD XF 366**



**345**



**350**



**400**



**355**



**358**



**361**



DISKOVI



Preporučeni broj obrtaja u minuti od 25.000 - 30.000 rpm

# TRI HAWK

FOR THE BUR CONNOISSEUR

Možete li da uklonite  
6 krunica za  
manje od 2 minuta?

10 + 1  
GRATIS

- ✔ Seče horizontalno i vertikalno
- ✔ Seče keramiku, metal, amalgam
- ✔ Utrošite manje vremena po pacijentu

PROBAJTE NAJBOLJE SEKAČE  
METALA I KERAMIKE NA SVETU




The logo for GHIMAS, featuring a stylized 'G' symbol followed by the brand name in a bold, sans-serif font.

**PRIVREMENI ISPUN  
VRHUNSKOG KVALITETA  
OD POZNATOG  
ITALIJANSKOG PROIZVOĐAČA**



**50g**  
PAKOVANJE



Dostupan u beloj i svetlo braon boji



MIX BORERA ZA BRUŠENJE

BORERI ZA KONZERVATIVU

POLUŽLEB

STEPENIK

(856) (850)  
(881) (882)

(837) (847k) (848)

(SG 835-010) (SG 836-014)  
(SG 830rl-010) (SG 830rl-012)  
(SG 830rl-014) (SG 801-010)  
(SG 801-012) (SG 801-018)  
(SG 801-021) (SG 801-025)





# Roto Dent

ROTIRAJUĆI DENTALNI INSTRUMENTI



060 30 93 080  
060 32 31 643

[www.rotodent.rs](http://www.rotodent.rs)  
[info@rotodent.rs](mailto:info@rotodent.rs)

[@roto\\_dent](https://www.instagram.com/roto_dent)

